

# Comfort Plus™ Cannulas

**Adult  
Pediatric  
Infant  
Neonatal  
Premature**



With the patient's ease of fit and maximum comfort in mind, the Comfort Plus line of Cannulas are designed lightweight and clear to ensure a comfortable fit for long-term use. The Over-the-Ear, "Bolo-Style" adjustable slide moves smoothly, yet provides a secure and comfortable fit while allowing maximum freedom of movement. The curved nasal prongs enhance and promote alignment of the cannula to fit naturally on the patient's face, eliminating the irritating flap required with conventional style cannulas. The Comfort Plus Soft Tip Line of Cannulas comes standard with "six-channel" safety tubing to ensure uninterrupted flow of oxygen to the patient. Also available in special lengths, please contact a Westmed representative for details.

Product Number	Description	Case Qty
<b>Comfort Plus Adult Soft Tip Cannula</b>		
0137	Micro w/7' Kink Resistant Tubing	.50
0194	4' Kink Resistant Tubing	.50
0197	7' Kink Resistant Tubing	.50
0198	14' Kink Resistant Tubing	.50
0199	14' Kink Resistant Tubing w/ Threaded Nut	.50
0195	25' Kink Resistant Tubing	.25
0196	50' Kink Resistant Tubing	.25
0190	Harness Only w/2" Oxygen Tubing Connector	.50
<b>Comfort Plus Pediatric Soft Tip Cannula</b>		
0183	7' Kink Resistant Tubing and Threaded Nut	.50
0187	7' Kink Resistant Tubing	.50
<b>Comfort Plus Infant Soft Tip Cannula</b>		
0173	7' Kink Resistant Tubing and Threaded Nut	.50
0177	7' Kink Resistant Tubing	.50
<b>Comfort Plus Neonatal Soft Tip Cannula</b>		
0163	7' Kink Resistant Tubing and Threaded Nut	.25
0162	7' Kink Resistant Tubing	.25
0167	7' Kink Resistant Tubing	.50
<b>Comfort Plus Premature Soft Tip Cannula</b>		
0152	7' Kink Resistant Tubing	.25
0153	7' Kink Resistant Tubing and Threaded Nut	.25
0157	7' Kink Resistant Tubing	.50

**EarMates™** are soft foam tubes designed to prevent ear chafing.

Product Number	Description	Case Qty
<b>Adult CO2 Cannula</b>		
0503	CO2/O2 Cannula w/male luer	.25
0504	CO2/O2 Cannula w/female luer	.25
<b>Demand Cannula</b>		
0507	Adult Demand Cannula w/4' Tubing	.25
<b>High Flow Cannula</b>		
0550	Adult High Flow Cannula w/7' Tubing	.25
0551	Adult High Flow Cannula w/25' Purple Tubing	.10
0552	Adult High Flow Cannula w/50' Purple Tubing	.10
<b>Cannula Mates</b>		
0553	Adult Soft Tip Cannula w/EarMates Attached	.25
<b>Cannula Accessories</b>		
0440	EarMates™ - charcoal, sold in pairs	50
0590	O2 Swivel Connector	25

Comfort Plus™ and EarMates™ are Trademarks of Westmed, Inc.

A safe and simple alternative for oxygen therapy patients.

## The BiFlo™ & PediFlo™ features and benefits:

- The BiFlo nasal mask is a simple and comfortable device designed as an alternative for oxygen therapy patients who cannot or will not tolerate a cannula or for those patients who have recurring necrosis.
- PaO2 and FIO2 values are comparable to those receiving oxygen with a conventional nasal cannula.
- The BiFlo delivers a gentle flow of oxygen. Unlike a nasal cannula, oxygen flow is not jetted into the nasal cavity.
- The PediFlo animated design reduces the apprehension generally encountered by younger patients.
- Latex Free.

**PediFlo Nasal Mask, Pediatric (Animated)**



**BiFlo Nasal Mask, Adult**

Product Number	Description	Case Qty
<b>BiFlo™ Adult Mask</b>		
0105	BiFlo Nasal Mask w/2" Oxygen Tubing Connector	.25
0101	BiFlo Nasal Mask w/7' Kink Resistant Tubing	.25
0101-25	BiFlo Nasal Mask w/25' Kink Resistant Tubing	.25
<b>PediFlo™ Pediatric Mask</b>		
0106	PediFlo Nasal Mask, Pediatric Animated w/7' Kink Resistant Tubing	.25

