


Engineered to meet the unique needs of your smallest patients



**Figure A**

**Figure B**

Westmed innovative solutions for NICU give the clinician the option to use the MiniHeart® Low-Flo with the Thayer Valved Tee (Figure A) or the VixOne™ with Infant Nebulizer Kit (Figure B).



**NICU Nebulizer features and benefits:**

- Westmed Solutions for Neonatal Intensive Care
- High-performance VixOne™ nebulizer with unique aerosol delivery system
- Very low dead volume (0.7cc)
- Nebulizes in upright positions or at angles up to 45 degrees
- The low-cost MiniHeart® Lo-Flo is ideal for flow sensitive ventilators
- The Neonatal Thayer Valved Tee™ reduces the amount of time that constant inhalation flow is interrupted
- The Neonatal Thayer Valve reduces cross-contamination, remains inline when not in use and less ventilator circuit manipulation

Part Number	Description	Case Quantity
<b>INFANT NEBULIZER KITS W/EXTENSION SETS</b>		
0236	VixOne™ Infant Nebulizer Kits .....	.50
100616	MiniHeart® Nebulizer w/Extension set .....	.25
<b>INFANT NEBULIZER KITS W/EXTENSION and THAYER VALVE</b>		
0234	VixOne™ Infant Extension w/Thayer Valve .....	.25
100611C	MiniHeart® Nebulizer w/Extension w/Thayer Valve .....	.25

Please log onto [www.westmedinc.com](http://www.westmedinc.com) for a complete listing of Westmed anesthesia and respiratory products.



Engineered to meet the unique needs of your smallest patients


Look to Westmed for innovative products that help you meet your **NICU** needs. From emergency ventilation, to aerosolized drug delivery to arterial and venous blood gas sampling, Westmed brings you a broad family of products to help maximize **patient care** and **comfort**, while helping you minimize costs.

Part Number	Description	Case Quantity
<b>1cc A line-syringe contains 25 units of balanced Heparin</b>		
3100-20	1cc .....	100
3100-25	1cc .....	100
<b>1cc A line with Blunt Tip Cannula-syringe contains 25 units of balanced Heparin</b>		
3100-30	1cc A line .....	100
<b>Cord Gas Kits-syringe contains 25 units of balanced Heparin</b>		
3182	2-1cc Syringes .....	100
<b>Comfort Plus Infant Soft Tip Cannula</b>		
0177	7' Kink Resistant Tubing .....	50
<b>Comfort Plus Neonatal Soft Tip Cannula</b>		
0167	7' Kink Resistant Tubing .....	50
<b>Comfort Plus Premature Soft Tip Cannula</b>		
0157	7' Kink Resistant Tubing .....	50
<b>Hyperinflation Sets</b>		
7-N7904-62	1/4L latex-free bag, 84" O <sub>2</sub> & pressure lines .....	12
7-N7905	1/2L latex bag, 84" O <sub>2</sub> & pressure lines .....	12
7-N7905-62	1/2L latex-free bag, 84" O <sub>2</sub> & pressure lines .....	12
7-N7901	1L latex bag, 84" O <sub>2</sub> & pressure lines .....	12
7-N7901-62	1L latex-free bag, 84" O <sub>2</sub> & pressure lines .....	12
7-N7910	Hyperinflation elbow, 84" O <sub>2</sub> & pressure lines .....	12
0499	Manometer (± 60 cm H <sub>2</sub> O) .....	1
0497	Manometer Cover .....	1
<b>700 Series- Ultra thin rotational molded</b>		
7-9704	Infant with inflation valve .....	20
7-9705	Neonate with inflation valve .....	20
<b>800 Series - Durable blow molded cushion</b>		
7-9804	Infant .....	20
7-9805	Neonate .....	20
<b>810 Series - Durable blow molded cushion with inflation valve</b>		
7-9814	Infant with inflation valve .....	20
7-9815	Neonate with inflation valve .....	20
<b>900 Series - Durable blow molded blue cushion</b>		
7-9904	Infant .....	20
7-9905	Neonatal .....	20

All masks are individually bagged



• A line and Cord gas kits



• Nasal cannula for babies: premature, neonatal and infants



• Hyperinflation sets



• Westmed disposable face mask

# Westmed