

**Aerosolized Medication Device for Infants. An alternative to aerosol "blow-by"**

**The PediNeb™** was designed to enhance the delivery of aerosolized medication, while reducing the apprehension encountered by young patients.

**PediNeb™ features and benefits:**

- Pacifier. Calms the infant/toddler during aerosol treatment. Medication is inhaled through the nose, not the mouth.
- Medication Port. Directs aerosolized medication to the nose and eliminates medication waste associated with the "blow-by" technique.
- VixOne small volume nebulizer. Clinically proven, providing superior medication delivery performance.
- Easy to operate, clean and efficient.
- Developed by a pediatric nurse practitioner.
- Latex free.



*The PediNeb™, a new modality in aerosolized medication delivery.*

Patent Number 5,904,140  
PediNeb™ is a Trademark of Westmed, Inc.

Product Number	Description	Qty
<b>PediNeb™ Nebulizer</b>		
0382	VixOne Nebulizer, Pacifier, 45 Degree Elbow, ..... 12 & 7' Kink Resistant Tubing	
0383	VixOne Nebulizer, Pacifier, Pediatric Tee, ..... 12 Poppel Extension Set & 7' Kink Resistant Tubing	
0385	Pedi-Neb, Pacifier with 45 Degree Elbow ..... 25	

