

Wound Irrigation Drainage Basin



The **SmartBasin** is exactly as it states... a SmartBasin. It is large enough to handle a sizable laceration on the leg or thigh yet easily positions to accommodate smaller wounds on the arm, hand or head. For quick and easy clean-up, the SmartBasin uses a removable disposable liner.

Product	Number	Qty
SmartBasin , Reusable	0575	10
SmartBasin Liners , Disposable	0576	25
SmartBasin 1 Reusable & 24 Disposables	0577	25

WinShield™ and SmartBasin™ are Trademarks of Westmed, Inc.

**To Control Exposure
from Fluid Pathogens
and Minimize
Post-Care Cleanup.**

Smart Basin features and benefits:

- **Lightweight.** Easy to handle, carry and store.
- **Disposable Liner.** Eliminates messy cleanups, laundry and sterilization costs.
- **Molded Hanger.** Does not take up additional floor space in the ER and allows for easy storage and access.
- **Angled.** Tapered basin directs all fluids to the drain spout and into a kick bucket.
- **Flexible Basin Positioning.** Indented on three of the four sides, the SmartBasin can be adjusted to accept any appendage regardless of patient-gurney position.

For your complete wound care needs, use the SmartBasin™ in conjunction with the WinShield™ wound irrigation products.
(see reverse side)

Wound Irrigation System with Splash Shield

Wound Irrigation System with Splash Shield Protection and Optional Irrigation Pressure

WinShield™

features and benefits:

- **Transparent Splash Shield.** Operator is protected from splatter while maintaining a clear view of the laceration.
- **Elongated Oval Shape.** Covers all or most of the laceration and directs the splattered solution away from the clinician.
- **Increased Fluid Dynamics.** Irrigation pressures from 4 to 12 psi effectively decrease bacterial contamination in treated wounds.
- **Side Bleed Ports.** Irrigation splatter is minimized by positioning the splash shield as close as possible to the wound area - The side bleed ports allow excess fluid to drain.
- **Auto-Fill System.** Spring loaded syringe is connected to a double check-valve, with 6' tubing for large volume irrigation.

For your complete wound care needs, use the WinShield in conjunction with the SmartBasin wound irrigation products.

(see reverse side)

Westmed



The **WinShield** is a transparent shell that provides the caregiver with universal precautions from the splatter of bio-hazard fluids. This shield is conveniently packaged alone or with a syringe.

<u>Product</u>	<u>Number</u>	<u>Qty</u>
WinShield , Splash Shield Only	0570	50
WinShield with 35ml Syringe	0571	40
WinShield with 20ml Syringe and Auto Refill System	0572	10
WinShield with 12ml Syringe	0573	10

WinShield™ and SmartBasin™ are Trademarks of Westmed, Inc.