

WHEN ACCURACY AND PATIENT PROTECTION  
ARE YOUR TOP OBJECTIVES, THE NEW  
BLOCKAIDE™ PULMONARY FUNCTION FILTER  
MEETS THE STANDARD



**Elliptical Mouthpiece**



**Round Mouthpiece**

# BLOCKAIDE™ Pulmonary Function Filter

The BlockAide™ filter media consists of a special blend of polymers with a highly stable electrostatic charge on every individual fiber, which allows it to consistently achieve high efficiencies by using both the electrostatic charge and mechanical impaction to remove airborne bacterial and viral particles. While being highly efficient in bacterial and viral filtration efficiency, the filter media still maintains a low airflow resistance to ensure diagnostic accuracy.

## BlockAide Provides:

- Meets ATS recommendations of Standardization of Spirometry for Diagnostic Spirometry Equipment with in-line filters having air flow resistance of less than 1.5cm H<sub>2</sub>O/L/S.
- Effective filtration of bacterial and viral particulates
- A family of sizes to meet every connection application
- Low flow resistance to ensure diagnostic accuracy
- Low dead space to ensure no hidden artifacts

Viral Filtration Efficiency .....	99.9%
Bacterial Filtration Efficiency .....	99.9%
Flow Resistance .....	.46 cm/H <sub>2</sub> O/L/S@14L/S
Dead Space Round MP .....	54 ml
Dead Space Elliptical MP .....	51 ml

**6901-1E & 6901-1R Machine Side ID/OD 27.23 / 30.18 for:** INCHECK, INCHECK DIAL, MULTISPIRO, CREATIVEBIOMEDICS, VITALOGRAPH, SCHILLER, MICRODIRECT, SDI PULMOGARD II, MOST COSMED, MIR, P.K. MORGAN 232, FUTURE MED SPIROBAND, PB-600, CDX SPIRO 850, SPIROMETRICS FLOWMATE, FLOWMATE+, SMI I, VACUMED, MICRO, MICRO PLUS, MICRO GP-GL, MICRO LOOP, PONY

**6902-2E & 6902-2R Machine Side ID/OD 30.48 / 33.15 for:** SENSORMEDICS, 2200, 2400, 2450, 2800, VMAX, KEYSTONE I, KEYSTONE II, PAIL, CYBER-MEDICS MOOSE, GOULD 5004, VIAYSIS NEW

**6903-3E & 6903-3R Machine Side ID/OD 45.47 / 47.45 for:** KOKO, KOKO LEGEND, KEYSTONE, KOKO-TREK, CYBERMEDIC FULL GAS SYSTEMS, COLLINS 4G & 3G FULL GAS SYSTEMS

**6904-4E & 6904-4R Machine Side ID/OD 26.26 / 28.73 for:** COLLINS DC-2

**6905-5E & 6905-5R Machine Side ID/OD 29.52 / 32.77 for:** ONE FLOW, JAEGER

**6906-6E & 6906-6R Machine Side ID/OD 33.63 / 37.39 for:** MEDGRAPHICS

## Filter Selection Chart

### Elliptical Patient Mouthpiece

<u>Part Number</u>	<u>Machine Side I.D. / O.D.</u>	<u>Case Quantity</u>
6901-1E	27.23 / 30.18	50
6902-2E	30.48 / 33.15	50
6903-3E	45.47 / 47.45	50
6904-4E	26.26 / 28.73	50
6905-5E	29.52 / 32.77	50
6906-6E	33.63 / 37.39	50

## Filter Selection Chart

### Round Patient Mouthpiece

<u>Part Number</u>	<u>Machine Side I.D. / O.D.</u>	<u>Case Quantity</u>
6901-1R	27.23 / 30.18	50
6902-2R	30.48 / 33.15	50
6903-3R	45.47 / 47.45	50
6904-4R	26.26 / 28.73	50
6905-5R	29.52 / 32.77	50
6906-6R	33.63 / 37.39	50

\*Part numbers ending in "E" have elliptical mouthpieces

\*Part numbers ending in "R" have round mouthpieces

\*Noseclips are included in every package

Distributed by:

**Westmed Inc.** • 5580 S. Nogales Hwy. • Tucson, AZ 85706

**Toll Free 800-975-7987 • Phone 520-294-7987 • Fax 520-294-6061**

Please visit us at [www.westmedinc.com](http://www.westmedinc.com) for a complete listing of Westmed products.

