

CO₂ Easy™

CARBON DIOXIDE DETECTOR

Seeing is
Believing




Westmed
When Your Care is Critical

INTRODUCING CO₂ Easy™

CARBON DIOXIDE DETECTOR

Optimize Visibility, Clarity and Confidence in CO₂ Detection

Featuring the largest visual indicator available, the CO₂ Easy™ is ideally mounted up front where clinicians and first responders can effectively view critical carbon dioxide levels without looking away from the patient. The non-toxic, pH sensitive indicator uses familiar purple-to-yellow color conversion when CO₂ is exhaled.



- Built in bacterial filter protects both patient and clinician.
- CO₂ Easy for patients weighing 15kg (33lbs.) or above.
- Pedi CO₂ Easy for patients weighing 1 - 15kg (2.2 - 33lbs.).
- 3 year shelf life.
- Designed for fast, easy airway assessment.
- Verifies proper endotracheal tube placement.
- Helps prevent unwarranted intubations.
- Confirms cardiac output during CPR.
- Expansive 2" Diameter Visual Indicator makes rapid breath-to-breath response easy to see.
- Anti-fogging agent reduces susceptibility to moisture and humidity, while ensuring continuous lens clarity. (CO₂ Easy only)
- Complies with Standard of Care Guidelines from the American Heart Association, American College of Emergency Physicians, American Association for Respiratory Care and the American Society of Anesthesiologists.



Part Number	Description	Case Quantity
562134	CO ₂ Easy Carbon Dioxide Detector	10
562135	Pedi CO ₂ Easy Carbon Dioxide Detector	10
562136	BagEasy w/ Adult Mask and CO ₂ Easy	6
562233	BagEasy w/ Peep, Adult Mask and CO ₂ Easy	6

CO₂ Easy is a trademark and Bag Easy is a registered trademark of Westmed, Inc.
© 2006 Westmed, Inc. All rights reserved. PN 74558/Rev. 05

Distributed by:

Westmed
When Your Care is Critical

Westmed, Inc.
5580 S. Nogales Hwy.
Tucson, AZ 85706
Phone: 800-975-7987
Fax: 520-294-6061
www.westmedinc.com

