

Hauge™ Airway Conscious Sedation Kit



Way Beyond an Airway...

Westmed
When Your Care is Critical

Hauge™ Airway for Total Airway Management

CONSCIOUS SEDATION KIT

The Hauge Airway is designed for simplicity and ease of use for most conscious sedation cases. The anatomically correct profile provides an immediate patent airway without stimulating a gag reflex. The Hauge Airway provides enough width to keep the tongue anterior in the oropharynx to ensure airway patency. Multiple applications for the disposable airway include:

- Awake Intubations
- Outpatient procedures
- Surgical MAC cases
- Upper Endoscopy procedures
- Bronchoscopy procedures
- Emergency room care
- For use with Esophageal Dilators
- Sleep Apnea Patients

Benefits for the provider:

- Immediate, established airway
- Universal O₂ and ETCO₂ ports
- Capability to deliver high flow O₂ orally to adult mouth breathers
- Large orifice accommodates a variety of scopes, suctioning
- Pre-procedure placement, with patient comfort in mind
- Device design accommodates elastic strap to ensure airway placement and stabilization
- Access to the oropharynx
- Bite block protection

Distributed by:



Can accommodate:

- Suction catheter
- Fiber-optic scope with an endotracheal tube
- Upper Endoscopy scope
- TEE scope

Additional Uses:

- Nebulizer treatment – Corrugated tubing is connected to the large center opening
- Airway patency after general anesthesia

CSK includes:

- Hauge Airway
- CO₂/O₂ supply line
- Head strap to secure airway
- Tongue blade
- 2" Pigtail with female lure

Ordering Information:

Part Number	Description	Case Quantity
6586	Large Airway CSK	10/bx
6587	Medium Airway CSK	10/bx
6595	Hauge II Airway CSK	10/bx

Westmed, Inc.

5580 S. Nogales Hwy • Tucson, AZ 85706
Phone : 800/975-7987 • Fax : 520/294-6061
www.westmedinc.com • sales@westmedinc.com

