

A TOTAL SOLUTION FOR YOUR  
DRY BUBBLE HUMIDIFIER  
REQUIREMENTS ACROSS THE  
ENTIRE CLINICAL FLOW RANGE



**Dry Bubble Humidifier  
For Low Flow Applications**



**Dry Bubble Humidifier  
For High Flow Applications**

# DRY BUBBLE HUMIDIFIER

Westmed's Family of Dry Bubble Humidifiers provides effective and consistent results for your patients across the entire clinical flow range. Both the "Low Flow" and "High Flow" Dry Bubble Humidifiers incorporate unique features to address your everyday requirements.



## Dry Bubble Humidifier For Low Flow Applications

- Pressure Relief Valve
- Humidifier Jar Engraved with Minimum/Maximum Water Level Points
- Interfaces with Oxygen Concentrator, Liquid O2 Systems, O2 Tanks and Wall Outlet Applications
- Home Care, Hospital & Alternate Care Applications
- Recommended Flow Rates of up to 6 LPM



## Dry Bubble Humidifier For High Flow Applications

- Pressure Relief Valve
- Humidifier Jar Engraved with Minimum/Maximum Water Level Points
- Unique Baffling System to Minimize Splash Exiting into the Tube
- Interfaces with Oxygen Concentrator, Liquid O2 Systems, O2 Tanks and Wall Outlet Applications
- Home Care, Hospital & Alternate Care Applications
- Recommended Flow Rates of up to 15 LPM

Distributed by:

<u>Description</u>	<u>Part Number</u>	<u>Case Quantity</u>
Dry Bubble Humidifier, Low Flow	0480	50
Humidifier Elbow Adapter Kit	0485	50
Dry Bubble Humidifier, High Flow	0795	25
Dry Bubble Humidifier, High Flow, Comfort Soft Plus™ Adult High Flow Cannula, w/7' Kink Resistant Green Oxygen Tubing	0792	10
Dry Bubble Humidifier, High Flow, Comfort Soft Plus™ Adult High Flow Cannula, w/14' Kink Resistant Green Oxygen Tubing	0793	10
Dry Bubble Humidifier, High Flow, Comfort Soft Plus™ Adult High Flow Cannula, w/25' Kink Resistant Green Oxygen Tubing	0794	10

The best of all Worlds! Westmed's Adult High Flow Comfort Soft Plus™ Cannulas are available for use with the Westmed High Flow Dry Bubble Humidifier. Maximum comfort for the patient while virtually eliminating ear ulceration combined with a High Flow Humidifier that vastly reduces splash and residual solution in your oxygen tubing. Flow rates are not compromised by restrictions and resistance caused by excess condensate in the tubing, a common occurrence with industry competition.



## High Flow Humidifier Proprietary Baffling System

### Adult Comfort Soft Plus™ Cannulas for High Flow Therapy

<u>Description</u>	<u>Part Number</u>	<u>Case Quantity</u>
Adult w/2" Connector & O2 Adapters	0548	25
Adult w/14' Kink Resistant Green Oxygen Tubing	0549	25
Adult w/7' Kink Resistant Green Oxygen Tubing	0550	25
Adult w/25' Kink Resistant Green Oxygen Tubing	0551	10
Adult w/4' Kink Resistant Green Oxygen Tubing	0554	25

**Westmed Inc.** • 5580 S. Nogales Hwy. • Tucson, AZ 85706  
**Toll Free** 800-975-7987 • **Phone** 520-294-7987 • **Fax** 520-294-6061  
 Please log onto [www.westmedinc.com](http://www.westmedinc.com) for a complete listing of Westmed anesthesia and respiratory products.

