

Uni-Lim™ & Pedi-Lim™

ADULT AND PEDIATRIC CO-AXIAL CIRCLE ANESTHESIA BREATHING CIRCUITS

The
Optimal
Choice



Westmed

When Your Care is Critical

Uni-Lim™ & Pedi-Lim™

ADULT AND PEDIATRIC CO-AXIAL CIRCLE ANESTHESIA BREATHING CIRCUITS

The Westmed Uni-Lim™ is the optimal choice for your Co-Axial Circle System needs. Westmed offers the right product at an affordable price.

- Designed to function as a standard circle anesthesia breathing system.
- Features a streamlined tube-within-a-tube design which eliminates twisting between the inspiratory and expiratory limbs.
- Eliminates bulky connectors at the patient interface.
- Provides natural heat and moisture conservation properties to that of a standard circle circuit.
- Engineered with a secure-fit inspiratory limb to prevent patient end disconnects. No floating internal tubes.
- ISO-standard connections for secure, leak-resistant attachments.
- Available in convenient Latex-Free Configurations.
- Available in custom configurations to meet your clinical needs.
- Complete anesthesia kits available.

Many Anesthesia Kit Options to Choose From

- Plain or ported elbows
- HME's, HMEF's, bacterial viral filters, and HEPA filters
- Standard and swivel gas sampling lines
- Internal and external gas sampling configurations
- Multiple anesthesia face mask options
- Guedel and Berman airways
- ECG leads
- Flex extenders with and without gas sampling
- Many other options available

Adult Breathing Circuits

Part Number	Description	Case Quantity
Uni-Lim 40" Adult Co-Axial Circle Anesthesia Breathing Circuit		

6656	GSE, 3L Latex Free Bag, Angled HMEF	20
------	-------------------------------------	----

Uni-Lim 60" Adult Co-Axial Circle Anesthesia Breathing Circuit		
---	--	--

6653	GSE, 1 Filter, 3L Latex Free Bag	20
------	----------------------------------	----

Uni-Lim 72" Adult Co-Axial Circle Anesthesia Breathing Circuit		
---	--	--

8041	GSE, 1 Filter, 2L Latex Free Bag	20
------	----------------------------------	----


Pediatric Breathing Circuits

Pedi-Lim 40" Pediatric Co-Axial Circle Anesthesia Breathing Circuit		
--	--	--

8121	GSE, 1 Filter, 1L Latex Free Bag	20
------	----------------------------------	----

Pedi-Lim 60" Pediatric Co-Axial Circle Anesthesia Breathing Circuit		
--	--	--

8122	GSE, 1 Filter, 1L Latex Free Bag	20
------	----------------------------------	----



Distributed by:

Ask your local representative about other standard and custom configurations that are available to meet your breathing circuit needs.

Westmed, Inc.

5580 S. Nogales Hwy • Tucson, AZ 85706
 Phone : 800/975-7987 • Fax : 520/294-6061
www.westmedinc.com • sales@westmedinc.com

